

Poster 4

A sequential treatment strategy for ex vivo profiling of live tumor fragments that mitigates tumor heterogeneity and tissue scarcity from core needle biopsies to characterize response to immunotherapies

Christina Scribano', Pichet Adstamongkonkul', Torey Browning', Sean Caenepeel', Nicholas Dana', Jackie Derrick', Thomas Dietz', Evan Flietner', Hilary Hernan', Christin Johnson', Nathan Marhefke', Payton McDonnell', Betsy Mulligan', Amreen Nasreen', Victoria Pope', Josh Porter', Abhijeet Prasad', Jordyn Richardson', Mikaela Schultz', Sidney Schneider', Mike Smith', Chetan Sood', Aishwarya Sunil', Lindsey Vedder', Ellen Wargowski', Hinco Gierman', Andreas Friedl', Laura Hrycyniak', Michael Korrer', T.S. Ramasubramanian'

¹Elephas, Madison, WI, ²University of Wisconsin-Madison, Madison, WI,



Introduction

- · Immune checkpoint inhibitors (ICIs) are highly effective treatment options for cancer, yet they fail to provide clinical benefit for most patients, primarily due to the low accuracy of approved biomarkers (eq. PD-L1 and Microsatelite Instability (MSI) / Mismatch Repair (MMR))
- Ex vivo cytokine profiling of live tumor samples has shown promise for improved prediction of response to αPD-1 blockade¹
- · This approach has been limited to tumor resections given the need for large amounts of tissue due to challenges associated with intra-tumor heterogeneity
- . Here, we present a novel sequential treatment strategy to detect T-cell response to ICI treatment in the limited tissue available from a single core needle biopsy (CNB)

Methods

Specimen processing: Surgical resections or CNBs were cut using specialized cutting instruments to obtain live tumor fragments (LTFs) which were dispensed into assay plates (resections = 200, 300 µm cubed fragments/well; CNBs (12-18 gauge) = ~12 fragments/well) and

Cross-well comparison of human tumor CNBs and surgical cross-well comparisor or numin furnor varies and surgical resections: LTFs were treated for 24 hours with either IgG or ICI (a PD-1). Conditioned media was then collected and cytokine concentrations were measured using a 30-plex kit from R&D Systems and run on a Luminex instrument.

Validation of the sequential treatment strategy: LTFs were treated with one of three regimens: (1) IgS for 28 hours, (2) acD3/acD28 for 28 hours or (3) IgG for 20 hours sidiowed by acD3/acD28 for an additional 8 hours. Conditioned media was sampled at 4, 20 and 28 hours for or yorkine profiles.


Enrollment of specimens from clinical trials: A total of 167 specimens were received for ex vivo profiling, 65 specimens were excluded due to shipping issues (n=4), sample processing issue (n=4), or having a form factor other than CNB (n=57) (eg, forceps biopsy or fine needle aspirate).

Ex vivo profiling of human tumor CNBs: Human tumor CNBs were Ex vivo profiling of human tumor CNBs: Human tumor CNBs were fragmented using a specialized outling instrument to create LTFs which were encapsulated in a proprietary hydrogel for ex vivo culture *LTFs comprising 145 samples from 102 human CNBs were treated for 20 hours with 1gG followed by ICI (dPD-1, n=70; aPD-1, n=6; aPD-1+a CTLA-4, n=26) for an additional 28 hours. Conditioned media was collected at -4, -20 and -48 hours for cytokine profiling.

Derivation of modified Z-scores and hierarchical clustering: A frimmed sample method was used to establish non-response to ICI transmet as a baseline usering he population median and medians because the sample of the production relationship of the sample measurements winch railed a U.C. brock were replaced by the popular median (i.e. 0). Values were also clipped to enforce a modified Z-score range of -20 to 20. Hierarchical clustering was then performed on the resulting sample (columns) and cytokine (rows) values. Ward's linkage method was used to determine cluster formation for both rows and


PD-L1 / MMR / MSI status: PD-L1, MMR, and MSI statuses were obtained from patient medical records. Where necessary, PD-L1 and MMR protein expression levels were characterized from treated histological sections collected at the end of the study period. Expression was then quantified by a board certified-pathologis

1 CNBs are incompatible with cross-well comparisons due to intra-tumor heterogeneity and limited available



A Schematic showing the distribution of tissue from a CNB in an experiment using cross-well configuration. B The percentage of cytokines out of a 30-cytokine panel with concentrations a lower limit of quantitation (LLOQ) for the cross-well configuration for resection LTFs (resection cross-well) and biopsy LTFs (CNB cross-well) following IgG control for 24 hours. There is no significant difference in the ability to detect cytokines between the two form factor protocols (p = 0.175). C Five cytokines from resection LTFs and biopsy LTFs plated in the cross-well configuration and treated with ICI for 24 hours. Significantly larger variability in cytokine concentrations between replicate wells of treated groups in biopsy LTFs compared to resection LTFs. ns, not significant, "p.c0.01; Tp.c0.0001.

2 A sequential treatment strategy that compares cytokine production between treated and control conditions within the same tissue addresses tumor heterogeneity




A Protocol to validate the sequential treatment strategy. Resection LTrs were used due to the availability of sufficient issue to address inter-web heterogeneity and enable meaningful cross-well comparisons. LTrs encapsulated in hydrogel were treated with corricle [Ga antiboty for 28 hours (gening, Occ390cCD28 for 26 hours to semitant 5 receils (bulby, or [56 for the first 20 hours and or CD340CD28 for the final 8 hours of culture (rec.) B Cyboline production rates (slope) following stimulation at 12 ours and an 120 were compared. C [First Audicion in 3 human trans theat and need [First S02] colored [GR339], uterine [UT3540]) was observed following oCD340CD28 stimulation (Stim) at 10 (bulb and 120 (red.) D Fold change in slope from Stim at 120 relative to slope from all 120 relation of the same stimulation (Stim) at 10 (bulb and 120 (red.) D Fold change in slope from Stim at 120 relative to slope from a similar and of induction for multiple cybidines (state of 13 decident capitalist induction rates) for 7 human tumors.

3 A sequential treatment strategy enables characterization of ICI response in CNBs where tissue is limited



Encapsulated LTFs from CNBs are sequentially treated with IgG control followed by ICI. Conditioned media is collected at the 4-, 20- and 48-hours to determine cytokine production rate during the control (4-20 hours) and ICI (20-48 hours) treatment phases. Fold change in the rate of cytokine production (slope of ICI treatment phase / slope of control treatment phase) is used

4 Differential cytokine response observed across 12 non-small cell lung cancer (NSCLC) CNB specimens



A LU4824 is a NSCLC CNB exhibiting an increased rate of Tiley and CXCL1 production following CIt teatment suggesting a potential cytokine response (multiple CNBs were obtained from this partner exacting profiting of registers with SLU4857 is a NSCLC CNB exhibiting a smaller rate of partner exacting profiting of registers with SLU4857 is a NSCLC CNB exhibiting a smaller rate response is unlikely (multiple CNBs were obtained from this patient enabling profiting of 2 registers which is a continuous control of the control of

Conclusions

- · A sequential treatment strategy for the assessment of response to immunotherapy is an effective approach to capture treatment-induced changes in cytokine response
- · Sequential treatment of CNBs circumvents issues of tumor heterogeneity inherent to cross-well comparison in settings of limited tissue availability
- · Unsupervised hierarchical clustering of cytokine profiling data collected with the sequential treatment strategy identifies a subset of specimens enriched for cytokine response following ICI treatment
- · This platform provides a scalable approach with the potential to change clinical practice for cancer patients being considered for treatment with immunotherapy

References

1. Voabil P, et al. Nat Med. 2021;27(7):1250-1261. doi:10.1038/s41591-021-01398-3 Iglewicz B, Hoaglin D. Vol 16. In: Mykytka EF, ed. The ASQC Basic References in Quality Control: Statistical Techniques. ASQC Quality Press; 1993.

3. Adstamongkonkul P, et al. AACR Annual Meeting 2025, Poster 6539, Session

Acknowledgments

- Some tissue samples were provided by the Cooperative Human Tissue Network, which is funded by the National Cancer Institute.
- The authors thank the University of Wisconsin Carbone Cancer Center BioBank, supported by P30 CA014500, for use of its facilities and services, and Mercy Research located at Mercy Sindelar Cancer Center led by PI Peter Di
- The authors admonstered by it is rest or the beautiful from the investig of Kansas Cancer Center's Biospecimen Repository Core Facility staff, funded in part by the NCI Cancer Center Support Grant P30 CA18842. The University of Louisville, Division of Surgical Cincology Biobank and The University of Tennessee Medical Center's Biobank for its service in providing biospecimens under the direction of Arigi Burubbank. MQ. PhD.
- The authors are grateful to the patients and participating sites in ELEPHAS-01 (NCT05478538), ELEPHAS-02 (NCT05520099) and ELEPHAS-04 (NCT06349642) for their contributions to this work
- Some figures were created in Biorender.com.

Scan here for more information



5 Unsupervised hierarchical clustering of cytokine profiling data from 145 CNB samples collected with the sequential treatment strategy identifies a subset of samples enriched for cytokine upregulation following ICI treatment

