

Label-free dynamic optical coherence microscopy assesses cytotoxicity in response to immune checkpoint inhibitors in ex vivo live tumor tissue

Chao J. Liu¹, Jewel Tulod¹, Yuan-I (Nina) Chen¹, Eric Wait¹, Christin Johnson¹, Ellen Wargowski¹, Todd Bakken¹, Sean Caenepel¹, Pierre Cunin¹, Andrew Durette¹, Evan Flitner¹, Michael Korrer¹, Jonathan N Ouellette¹, Christina Scribano¹, Max Schiller¹, Lindsey Vedder¹, Kevin Eliceiri² and Michael Szulcowski¹

¹Elephas Biosciences, Madison, WI, ²University of Wisconsin-Madison, Madison, WI.



Poster 113

Introduction

- Cancer imaging assays typically use external labels, such as dyes, to assess treatment induced cytotoxicity
- Although dyes work well with tumor cell line monolayers and spheroids, they often fail to prove adequate for primary tumor tissues, where the native tumor microenvironment and architecture are maintained
- Dynamic optical coherence microscopy (dOCM), an emerging noninvasive, label-free imaging modality, analyzes temporal signals for fluctuations due to metabolic processes and cellular motions
- dOCM allows for longitudinal and nondestructive assessment of cell viability deep in intact tissue, removing the need for dyes
- Here we use dOCM to detect cytotoxic changes in tumor tissue in response to immune checkpoint inhibitor (ICI) treatment, and assess its alignment with the ground truth dye propidium iodide (PI)

Methods

Sample preparation: CT26 tumors were cut into fragments using an automated cutting instrument (Elephas). The tissue fragments were embedded in Matrigel (Corning) in polytetrafluorethylene standing well inserts (Millipore Sigma) which were positioned in multi-well glass bottom plates (Ibidi).

Imaging: A customized multimodal imaging system (Elephas) was used to obtain dOCM and multiphoton microscopy (MPM) imaging data with modalities coaligned through a 1.8x, 0.8 numerical aperture objective.

dOCM imaging of cytotoxicity: Images were acquired using Elephas custom acquisition software. Time resolved, volumetric dOCM data were obtained to extract structural and viability data. To observe ICI-mediated cytotoxicity, fragments were imaged every 12 hours, except for during the addition of treatment. Baseline data were collected during the first 12 hours (-24 to -12 hours of the experiment). Treatments (IgG, aPD-1 or aPD-1 + aLAG3) were added at the 0 hour time point.

MPM imaging of cytotoxicity: A titanium:sapphire laser (Coherent) was tuned to 740 nm to excite nicotinamide adenine dinucleotide (phosphate) (NAD(P)H) and PI (1 µg/mL, Invitrogen). PI was added to the culture media 4 hours prior to imaging. The fragments were imaged at 0, 24 and 48 hours after treatment.

In vivo tumor growth inhibition: BALB/c mice implanted with CT26 tumors were treated with PBS (vehicle control) or 100 µg aPD-1 or aPD-1 + aLAG3 (N=30 per group) starting five days post implantation and continuing twice every week until tumors in the control group reached 2000 mm³.


Immunohistochemistry: Immunohistochemistry was performed on representative CT26 tumors from each of the treatment groups harvested 26 days after implant. Tumors were formalin-fixed, paraffin-embedded (FFPE), sectioned and stained with CD8, cleaved caspase 3 (CC3) and Ki67 antibodies.

Immunofluorescence: Multiplex immunofluorescence was performed on representative CT26 tumors from each of the treatment groups at the end of experimentation. Tumors were formalin-fixed, paraffin-embedded (FFPE), sectioned and stained with CD8 and CD4 antibodies, and DAPI.


Ex vivo cytokine analysis: Untreated, subcutaneous CT26 tumors were harvested at 200-300 mm³. Core needle biopsies (CNB) were collected from excised tumors and cut into 300 µm thick slices. Slices were treated with IgG, aPD-1 or aPD-1 + aLAG3 antibodies for 20 hours. Conditioned media was collected for IFNγ analysis by ELISA after 20 hours of treatment.

1 dOCM detects chemically-induced cytotoxicity in syngeneic tumor tissue


Use of dOCM to assess cell viability in CT26 tumor fragments



Accuracy of dOCM in assessing H₂O₂-induced cytotoxicity in CT26 tumor fragments using PI staining as ground truth





2 dOCM detects heat-shock-induced cytotoxicity in primary human bladder tumor tissue




dOCM imaging of a 300 µm thick primary human bladder tumor fragment heat-shock treated for 12 hours. (A) Images were collected prior to treatment (left) and after 10 hours of treatment (right). Scale bar = 100 µm. A qualitative reduction in dOCM contrast can be seen following treatment which is represented quantitatively in (B) as a change in the percentage of live voxels (live voxels/structural voxels) over time relative to time point 0 hrs. A timeline video of the experiment and a supplemental experiment performed in primary human lung tumor tissue may be viewed by scanning the QR code on the bottom right hand corner of this poster.

3 dOCM detects ICI-mediated cytotoxicity in CT26 tumor fragments ex vivo




T-cells remain within fragments and IFNγ is increased with ICI treatment ex vivo




Time course dOCM imaging analysis of CT26 tumor fragments treated with immune checkpoint inhibitors. (A) Schematic showing experimental protocol of tumor processing, treatment and imaging. (B) Representative dOCM images of CT26 tumor fragments prior to treatment and following 48 hours of treatment with aPD-1 or IgG. Reduced dOCM contrast can be seen in the aPD-1-treated fragment at 48 hrs after treatment. Scale bar = 100 µm. (C) Time course dOCM imaging analysis of CT26 tumor fragments treated with IgG, aPD-1 or aPD-1 + aLAG3 (N = 4 tumors with a range of 12 fragments per time point and treatment). Relative live voxels (%) was calculated by dividing the number of live voxels by the total number of structural voxels over time relative to time point -24 hrs. A two-way ANOVA revealed a significant interaction between day and treatment (F(10,18) = 3.430, P = 0.0006). Tukey's multiple comparisons analysis showed a significant difference in the aPD-1 treated group between -24 hrs and 48 hrs (P=0.009). Error bars represent SEM. *P<0.05, **P<0.01, ***P<0.0001. (D) Representative multiplex immunofluorescence images for CD8 (orange) and CD4 (green) staining show T-cells are maintained in fragments from all treatment groups. Scale bar = 50 µm. (E) Graph showing the fold change from IgG in IFNγ cytokine concentrations for aPD-1 (100 µg/µL) or aPD-1 (50 µg/µL) + aLAG3 (50 µg/µL) following 20 hrs of treatment. N = 6 and 7 tumors for aPD-1 and aPD-1 + aLAG3 groups, respectively. Error bars represent SEM.

4 ICI-mediated CT26 tumor growth inhibition in vivo aligns with dOCM observations ex vivo



(A) In vivo CT26 tumor growth curves from mice treated with vehicle, aPD-1, or aPD-1 + aLAG3 (N = 30 per group). Error bars represent SEM. Significant tumor growth inhibition was observed with aPD-1 alone or the combination of aPD-1 + aLAG3. (B) Immunohistochemistry conducted on tumors harvested 26 days after implant showing qualitative increases in CD8 and CC3 staining following treatment with aPD-1 or aPD-1 + aLAG3. Scale bar = 50 µm.

5 MPM imaging of PI detects ICI-mediated cytotoxicity in CT26 tumor fragments ex vivo



(A) Representative MPM images showing IgG or aPD-1 + aLAG3 treated tumor fragments at 24 hrs of treatment compared to 48 hrs of treatment. (B) Bar graphs report the difference in dead cell count/mm² at 48 hours of treatment versus 24 hours of treatment. The difference in dead cell count/mm² was calculated by dividing the number of PI-positive cells by the NAD(P)H total area. n = 6 fragments per group across 2 tumors with 4 z-planes imaged per fragment. Error bars represent SEM.

Conclusions

- dOCM, a non-destructive, label-free imaging technique, can be used to measure ICI-mediated cytotoxicity in tumor tissue
- Ex vivo assessment of ICI-mediated cytotoxicity with dOCM reveals efficacy that is aligned with in vivo tumor growth inhibition in the murine syngeneic tumor model CT26
- dOCM and the ground truth dye propidium iodide correspond in their assessment of ICI-mediated cytotoxicity ex vivo
- These data support the use of dOCM to assess the efficacy of immunotherapy in tumor tissue ex vivo

Acknowledgments

Human tissue samples were provided by the Cooperative Human Tissue Network which is funded by the National Cancer Institute. The authors thank the following individuals for their contributions to this work: John R Baffert, Nicholas Kapilewski and Katherine Hexroad for their advice in the development of this poster. Some figures created in BioRender.com.

References

- Riccardi C, Nicotri L. *Nat Protoc*. 2006;1(3):624-640. doi:10.1038/nprot.2006.238.
- Rosenberg M, Azavedo JF, Isler A, et al. *Biomol Opt Express*. 2024;15(1):1408-1417. Published 2024 Feb 5. doi:10.1364/boe.2023.1501408.
- Liu CJ, Smith TJ, Wang Y, et al. *Biomol Opt Express*. 2024;15(1):1408-1417. Published 2024 Feb 5. doi:10.1364/boe.2023.1501408.
- Amato SM, An AS, Lachin JR, et al. *Cancer Discov*. 2024;14(2):228-240. doi:10.1158/2156-8420.CCR23-0833.
- Song AT, Borczak A, Smith D, et al. *NPJ Precis Oncol*. 2023;7(1):127. Published 2023 Dec 8. doi:10.1038/s41601-023-00074-3.
- Shimokawa M, Kishikawa T, Rish S, et al. *Cancer Sci*. 2020;111(4):1228-1240. doi:10.1111/cas.14320.
- Benayahu A, Elomayor L, Lior H, Kivimäki M, Kivimäki M. *Cancer Sci*. 2023;114(3):715. Published 2023 Mar 13. doi:10.1111/cas.14320.
- Gilman P, Vasta V, Kivimäki M, Blum M, Blum M. *Assay Drug Dev Technol*. 2015;13(5):529-546. doi:10.1089/adt.2014.0018.
- Saito C, Yoshimura K, Kivimäki M, Kivimäki M. *Sci Rep*. 2023;13(1):2553. Published 2023 May 31. doi:10.1038/s41598-023-28271-2.
- Ortiz E, Gray SG, Finn SP. *Methods Mol Biol*. 2023;2465:265-275. doi:10.1007/978-1-0716-3056-5_17.

Scan here for more information

